

INVACARE®
GET-U-UP™
HYDRAULIC STAND-UP LIFT



Comfortable and secure, the **Invacare® Get-U-Up™ Stand-Up Lift** is the answer to your patient transfer needs. Its ergonomic styling and easy maneuverability make the Get-U-Up Lift an ideal product for everyday patient handling. The Get-U-Up Lift helps ensure caregiver security and patient peace of mind. Ideal for use with partial weight-bearing patients as well as those needing rehabilitation support, the Invacare Get-U-Up stand-up lift offers safety, comfort and stability.

Features

Innovative adjustment features adapt to a wide variety of body sizes and shapes. That means comfort at any height, width and weight range up to 350 lb.

Low friction casters significantly improve rollability and contribute to caregiver safety

Multi-function slings allow for quick toileting, stand-assist or full-support seated transfer

Adjustable leg support features a soft pad with lateral contour that secures legs in multiple positions for comfort at any height

Retractable, non-slip footplate provides stable base for transfer


Under-bed clearance of 4.5" allows accessibility to most homecare-style beds

Ergonomic features make it ideal for toileting transfers, as well as dependent and rehabilitation purposes

Fast, simple sling attachment that protects against accidental disengagement; wide push handle, good access for clothing removal; and comfortable slings with extended back support that prevent pressure in the under arm area

Pump handle can rotate from side to side for convenience of caregiver

HCPCS Code: E0630



1. Standing sling* Model # R130/R131

2. Transfer sling* Model # R134/R136

*Complete instructions for lift usage and sling procedures are found in the Invacare operating manual. Only Invacare® slings and accessories should be used on Invacare® lifts.



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2013
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-249
Rev. 12/13
131046



Model #	Height At Sling Hook Up	Base Dimensions	Casters	Weight	Weight Capacity	Low Battery Alarm	Emergency Controls
GHS350	36"-65.7" H	4.5"H Clearance 37"W Open 23"W Closed 40.5"L	3" Front 5" Rear	88 lb.	350 lb.	N/A	N/A



Battery	Charge Time	Charger Input	Charger Output	Limited Warranty
N/A	N/A	N/A	N/A	3 Years Lift 1 Year Pump



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.