PHILIPS

RESPIRONICS

DreamStation Go

Portable PAP therapy system





© 2017 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice.

www.philips.com/respironics

Caution: US federal law restricts these devices to sale by or on the order of a physician.

RRDPGH EB 3/1/17 PN 1133355 MCI 4107883 1010 Murry Ridge Lane, Murrysville, PA 15668 800-345-6443 • 724-387-4000

Convenience without compromise





Micro-flexible tubing

The **DreamStation Go** is a compact and highly portable PAP system designed for patients who refuse to compromise.





Powerful, portable, performance

DreamMapper cross-device connectivity

Users of DreamMapper, our patient self-management tool, ca sync data from multiple DreamMapper-enabled devices* to one account. Users have the ability to track and monitor their combined progress across both home and portable PAP devi-in one convenient dashboard on their smartphone, tablet, or computer

Micro-rtexible tubing DreamStation Go features the smallest PAP tubing we've ever of Engineered to minimize bulk and improve mobility, our 12mm tul is 58% smaller than our 22mm performance tubing. It offers enhanced flexibility, portability, and connectivity. • The click-connect PAP-side cuff simplifies connection of the unself the tet PAP.

Innovation for people on the go

DreamStation Go provides the portability and packability travelers want and the effective therapy they need.

65% smaller

The DreamStation Go is 65% smaller than ResMed's AirSense 10 and AirStart 10.

48% lighter*

The DreamStation Go is 48% lighter than ResMed's AirSense 10 and AirStart 10.

33% fewer components*

The DreamStation Go requires 33% fewer components to pack than its closest competitors.

- *2016 internal measurement of ResMed's AirSense 10 and AirStart 10 with side cover installed. Includes power supply Data on file.
- **When compared to HDM Z1 and Somnetics Transcend devices, (2016).

Powerfully convenient features



Onboard USB port allows for easy charging of mobile devices.

*Sold separately.

**RTCA/DO-160G section 21, category M.





Bluetooth is standard with every device.



Complies with FAA requirements for in-flight use.**



Device is compatible with our DreamMapper patient self-management tool.

12mm tube

Our new-to-market 12mm micro-flexible tubing is our smallest and lightest ever.

Click-connect

Click-connections on each end of the micro-flexible tubing make connecting and disconnecting the device and mask simple.

Accessories create a world of opportunities



Overnight battery Compact and lightweight with a recharge time of just 4-5 hours.



Small travel kit

Medium travel kit Kit includes an easy-to-Kit includes a protective pack protective case for case that accommodates the DreamStation Go and the DreamStation Go a separate, washable bag and overnight battery. the mask and tube. Also includes separate, washable bag for the

k and tube



Air filters Your choice of a reusable pollen or disposable 30-day fine filter.



DreamStation Go Auto specifications

	4 to 20 cm H_2^0
Pressure mode(s)	CPAP, APAP
Flex pressure relief	1 to 3
Flex modes	
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable, Smart Ramp
Dimensions	15 cm L x 15 cm W x 5.8 cm H
	844 g (1.86 lb.)
	Reusable pollen, combination fine
Filters Device setup	Reusable pollen, combination fine Touchscreen color display
Device setup	Touchscreen color display On-board: > 1 year
Device setup Data storage capacity	Touchscreen color display On-board: > 1 year MicroSD card (optional): lifetime
Device setup Data storage capacity Compliance meter	Touchscreen color display On-board: > 1 year MicroSD card (optional): lifetime Breathing detection

Ordering information

Flow Generators	Part number	
DreamStation Go CPAP	DSG400S11	
DreamStation Go Auto CPAP	DSG500S11	
Power	Part number	
DreamStation Go overnight battery	1133281	
Travel kits	Part number	
DreamStation Go small travel kit*	1133275	
DreamStation Go medium travel kit**	1133276	
Consumables/Replacement	Part number	
DreamStation micro-flexible 12mm tubing		
Power cord 6ft white	1132551	
Power cord 10ft white		
Beauty panel - front		
Beauty panel - rear	1133742	
Reusable pollen filter 2-pk		

DreamStation Go CPAP specifications

Pressure range	4 to 20 cm H ₂ O
Pressure mode(s)	CPAP
Flex pressure relief	
Flex modes	
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable, Smart Ramp
Dimensions	15 cm L x 15 cm W x 5.8 cm H
Weight	844 g (1.86 lb.)
Filters	Reusable pollen, combination fine
Filters Device setup	Reusable pollen, combination fine Touchscreen color display
Device setup	Touchscreen color display On-board: > 1 year
Device setup Data storage capacity	Touchscreen color display On-board: > 1 year MicroSD card (optional): lifetime
Device setup Data storage capacity Compliance meter	Touchscreen color display On-board: > 1 year MicroSD card (optional): lifetime Breathing detection

*Small travel kit includes an easy-to-pack protective case for the DreamStation Go and a separate washable mask and tube bag.

"Medium travel kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes a separate washable mask and tube bag.

Learn more

For more information about DreamStation Go, visit www.philips.com/dreamstationgo.

