Hill-Rom.

NP50 Prevention Surface



An economical prevention mattress for daily patient needs

The Hill-Rom non-powered NP50 prevention surface is a high quality and affordable foam mattress that addresses pressure redistribution. This surface, designed for use in acute care and extended care environments, is constructed of high quality, medical grade foam. This lightweight mattress is also durable and reversible from head to foot.

At a glance

- Single density foam prevention mattress
- Fire barrier included in all surfaces
- Lightweight, easy-to-clean, wipeable cover
- Fluid-resistant, leak proof fabric protects the surface from fluid penetration
- Durable ticking accommodates heavy use in a rigorous environment
- Anti-skid bottom helps keep surface in place
- Enclosed design with no protruding parts
- EPA-registered anti-microbial agent in foam to provide additional protection from odor-causing bacteria
- Optional bolster/raised foam perimeter
- Reversible from head to foot
- Articulation cuts allow the surface to conform to the bed frame
- Accommodates weight up to 350 lbs*
- * Bariatric versions available, accommodating up to 700 lbs
- * Psychiatric versions available



Hill-Rom_®

NP50 Prevention Surface



Hill-Rom_® NP50 Prevention Surface Specifications – Standard Mattress

Single layer foam construction


Firm support; 1.8 PCF
55 ILD; Density 1.8 PCF
78", 80" and 84"
35.5" +/-3.5"
6"
Approx. 24 lbs
350 lbs
29, CPSC 16CFR1633, Boston IX-11
2 years full
300 gm/sqm/24 hrs.

Hill-Rom_® NP50 Prevention Surface Specifications – Bariatric Mattress

Single layer foam construction

Base layer: 7"	Firm support; 1.8 PCF
Firm support base	46 ILD; Density 1.8 PCF
Length	80" and 84"
	39", 48"
Height	7"
Weight	Approx. 35 lbs
Weight limit	700 lbs
Flammability	Cal129, CPSC 16CFR1633, Boston IX-11
Warranty	2 years full
MVP	300 gm/sqm/24 hrs.


NP50 Flat deck surface




Head of Bed 30° 116 lb female



Head of Bed 30° 148 lb female



Head of Bed 30° 315 lb male



Representative sample for supine patient at weight indicated. Pressure results will vary depending on patient weight and position. Always use in conjunction with facility protocols for patient repositioning and pressure ulcer prevention.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

