

QUICKIE[®]

Power Portfolio





Table of Contents

- Zippie® ZM-3102
- Quickie® QM-7 Series3
- Quickie® Pulse Series6
- Quickie® S-6 Series8
- Quickie® P-222 SE10
- PG Drives Technology12
- Back Upholstery Options12
- Electronics13
- Base Options16
- Accessories17
- Armrests18
- Positioning Options18
- Legrests19
- Built-4-Me20
- JAY® Seating Selection20

Why Quickie?

You drive our design. Quickie wheelchairs are designed for you and by you. In development, we consider how you interact with your wheelchair to maintain maximum efficiency, performance, fit, and style. This results in the creation of products that are innovative, durable, versatile, and customizable. That is the Quickie Advantage.

Innovation. We are continually developing new products, improving our flagship products, and offering new options and solutions to enhance your lifestyle.

Quality and durability. Quickie wheelchairs undergo rigorous testing that surpasses industry standards, so you can be confident that your product will hold up and perform throughout its lifetime.

Versatility. We offer a comprehensive portfolio of products to meet a wide range of individual lifestyles and clinical needs. Once you choose a Quickie, our goal is that it continues to be the right chair for you. We offer maximum interchangeable parts and adjustability in our designs to help accommodate changes in condition, environment, or lifestyle that you may experience.

Built-4-Me. Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint.

zippie® ZM-310™

The Zippie ZM-310 combines our smallest, most reliable power wheelchair base with clinically proven pediatric seating and positioning options to provide children maximum accessibility and independence.



Compact Base with Serious Power!

- 24" wide, 35" long base for indoor accessibility
- 19" turning radius for intuitive driving and maneuverability
- 13" drive wheels for climbing obstacles and gliding over terrain

Positioning Options

- 50° of power tilt
- 50° of manual tilt
- 40° of manual recline
- 9" or 12" of power seat elevate



Pediatric Seating

- Available in sizes ranging from 10" x 10" to 16" x 19"
- ASAP II seat frame offers in-field width, depth, armrest, legrest, and back height adjustability
- Angle-adjustable centermount footrest with pediatric footplate and 3" to 13" lower leg lengths

ZIPPIE ZM-310 SPECIFICATIONS

- HCPCS Code**
K0848, K0856
- Drive Wheel**
Mid Wheel
- Max Speed**
6 mph
- User Weight Capacity**
165 lbs
- Width**
10" to 16"
- Depth**
10" to 19"
- Seat Height**
15.5" to 21.8"
- Overall Base Width**
24"
- Overall Length**
35"
- Turning Radius**
19"
- Drive Wheels**
13" with flat-free insert
- Transit Tested**
WC-19

Live Without Limits

The Quickie QM-7 series offers a high-performing, reliable, and stylish mid-wheel drive power solution for a client who qualifies for a Group 3 Standard (QM-710), Group 3 Heavy Duty (QM-715 HD) or Group 4 High Speed (QM-720) power mobility base.

QUICKIE® QM-7 Series



Available in Nine Group 3 and Group 4 Models:

GROUP 3: Standard Wt. Cap.

- Quickie QM-710 BC (K0848):
No Power Option Base with ASAP II Rehab Seat
- Quickie QM-710 SC (K0856):
Single Power Option Tilt (SPOT) or Recline Base with ASAP II or Recline Rehab Seat
- Quickie QM-710 MPC (K0861):
Multiple Power Option Base (SPOT Tilt and Recline) with 3.7 Recline Seat

GROUP 3: Heavy Duty

- Quickie QM-715 HD BC (K0850):
No Power Option Base with ASAP II Rehab Seat
- Quickie QM-715 HD SC (K0858):
Single Power Option Tilt (SPOT) with ASAP II Rehab Seat
- Quickie QM-715 HD MPC (K0862):
Multiple Power Option Base (SPOT Tilt and Recline) with 3.7 Recline Seat

GROUP 4: High Activity

- Quickie QM-720 BD (K0868):
No Power Option Base with ASAP II Rehab Seat
- Quickie QM-720 SD (K0877):
Single Power Option Tilt (SPOT) or Recline Base with ASAP II or Recline Rehab Seat
- Quickie QM-720 MPD (K0884):
Multiple Power Option Base (SPOT Tilt and Recline) with 3.7 Recline Seat

QUICKIE® QM-7 Series



Power Base Performance

- Linix 5 mph, 6 mph, or 8 mph, 4-pole motor packages
- Intuitive mid-wheel drive handling
- Contoured, sleek base available in 16 shroud color options
- SpiderTrac™ Suspension for superior traction and handling across multiple terrains
 - Gas dampening system for smooth transitions
 - 3" curb climbing capability
 - 22° transition angle helps prevent high centering on curb cuts



Seating Options

ASAP II Seat Frame

- In-field width, depth, armrest, legrest, and back height adjustability
- Wide range of seat sizes
- Manual reclining backrest for increased function and sitting tolerance
- Optional fold down backrest for easy transport
- Extensive seating and positioning options


3.7 Recline Seat Frame

- 14" – 22" seat widths and depths
- Ultra low seat-to-floor heights – as low as 16" when combined with tilt and elevate!
- Integrated JAY J3 back for seamless functionality with power recline and shear reduction system
- Cantilever Reclining Armrests for optimal positioning at any back angle
- Low profile electronics enclosure shell for a clean look and compatibility with accessories



Power Options

- NEW 3.7 Power Recline with Power Shear Reduction
- 50° of Power Tilt to optimize positioning
- 12" and 9" Power Seat Elevate for increased access to cabinets and interaction with peers
- Power Elevating Legrests with Articulation
 - Swing-Away
 - Centermount with Single/Dual Footplate for angle adjustability as a single footplate or separate footplates



NEW! Options for Consumer Convenience

- **RAM® Phone and Tablet Mounts** allow mobile phones and/or tablets to be connected and disconnected for easy viewing
- **Active Controls Power Buddy** provides an auxiliary USB power source for charging devices and optional LED lights for improved visibility
- **RAM® Mounts Self-Leveling Cup Holder** uses gimbal system to keep drinks level during travel or when stopped on inclines
- **Active Controls ReCharge™ with Socket** eliminates the challenges of connecting and disconnecting the battery charger by providing an easy-to-reach adapter

Quickie QM-7 Series Specifications

| SPECIFICATIONS | QM-710 BC | QM-710 SC | QM-710 MPC | QM-715 HD BC | QM-715 HD SC | QM-715 HD MPC | QM-720 BD | QM-720 SD | QM-720 MPD |
|-----------------------|-------------------|------------------------|-------------------|-------------------|------------------------------------|-------------------|-------------------|------------------------|-------------------|
| HCPCS Code | K0848 | K0856 | K0861 | K0850 | K0858 | K0862 | K0868 | K0877 | K0884 |
| Drive Wheel Position | Mid | Mid | Mid | Mid | Mid | Mid | Mid | Mid | Mid |
| Speed | up to 6.0 mph* | up to 6.0 mph* | up to 6.0 mph* | 5.0 mph | 5.0 mph | 5.0 mph | 8.0 mph | 8.0 mph | 8.0 mph |
| Overall Weight | 310 lbs. | 350 lbs. | 375 lbs. | 310 lbs. | 350 lbs. | 375 lbs. | 310 lbs. | 350 lbs. | 375 lbs. |
| Overall Width | 25" | 25" | 25" | 25" | 25" | 25" | 25" | 25" | 25" |
| Overall Length | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests | 36" w/o footrests |
| User Weight Capacity | 300 lbs. | 300 lbs. | 300 lbs. | 400 lbs. | Tilt: 400 lbs. Recline: 350 lbs | 350 lbs. | 300 lbs. | 300 lbs. | 300 lbs. |
| Seat Widths | 12" to 22" | 12" to 22" | 14" to 22" | 16" to 24" | 16" to 24" | 16" to 22" | 12" to 22" | 12" to 22" | 14" to 22" |
| Seat Depths | 12" to 22" | 12" to 22" | 16" to 22" | 16" to 22" | 16" to 22" | 16" to 22" | 12" to 22" | 12" to 22" | 16" to 22" |
| Seat To Floor Height | 16.3" to 21" | 16.1" to 20.6" | 16" to 19.3" | 16.3" to 20.3" | 16.1" to 20.8" | 16.6" to 19.6" | 16.3" to 21" | 16.2" to 20.6" | 16" to 19.3" |
| Seat Frame Options | ASAP II | ASAP II or 3.7 Recline | 3.7 Recline | ASAP II | ASAP II or 3.7 Recline | 3.7 Recline | ASAP II | ASAP II or 3.7 Recline | 3.7 Recline |
| Power Seating Options | N/A | Tilt or Recline | Tilt & Recline | No Power Options | Tilt or Recline | Tilt & Recline | No Power Options | Tilt or Recline | Tilt & Recline |

* Optional 5.0 mph and 8.0 mph motor packages also available

Live Without Limits

The Quickie Pulse is a compact, durable power wheelchair that offers a wide range of effective seating and electronics options to meet the needs of Group 3 users.

QUICKIE® Pulse™



Available in Six Models:

- Quickie Pulse 5 BC (K0848): 5 MPH, 2-pole motor, No Power Option base with ASAP II Rehab Seat
- Quickie Pulse 6 BC (K0848): 6 MPH, 4-pole motor, No Power Option base with ASAP II Rehab Seat
- Quickie Pulse 5 CC (K0849): 5 MPH, 2-pole motor, No Power Option base with Captain's Seat
- Quickie Pulse 6 CC (K0849): 6 MPH, 4-pole motor, No Power Option base with Captain's Seat
- Quickie Pulse 6 SC (K0856): 6 MPH, 4-pole motor, Single Power Option Tilt (SPOT) or Power Recline with ASAP II Rehab Seat or 3.7 Recline Seat
- Quickie Pulse 6 MPC (K0861): 6 MPH, 4-pole motor, Single Power Option Tilt (SPOT) and Power Recline with 3.7 Recline Seat



Compact Power Base

- Intuitive mid-wheel drive handling
- Small footprint and compact design for indoor accessibility
- Flexible, independent suspension system for moderate outdoor performance
- Linix 2- or 4-pole motors



Seating Options

Captain's Seat (CC)

- Convertible seat pan for flexibility in seating options
- Semi-recline and fold-down backrest
- Seat size, legrest, and armrest adjustability
- Mounting capabilities for O2, crutch, or walker holders

3.7 Recline Rehab Seat (SC & MPC)

- Option of Recline Curved back or Integrated JAY J3 Back for seamless functionality with power recline and shear reduction system
- Seat-to-floor heights as low as 16.3"
- Cantilever Reclining Armrests stay parallel to ground when reclined



ASAP II (BC & SC)

- In-field width, depth, armrest, legrest, and back height adjustability
- Manual reclining backrest for increased function and sitting tolerance
- Optional fold down backrest for easy transport

Power Options

- 3.7 Power Recline with Power Shear Reduction
- 50° of Single Power Option Tilt (SPOT) to optimize positioning
- 9" Seat Elevate for increased access to cabinets and interaction with peers
- Power Elevating Legrests with Articulation



NEW! Options for Consumer Convenience

- **RAM® Phone and Tablet Mounts** allow mobile phones and/or tablets to be connected and disconnected for easy viewing
- **Active Controls Power Buddy** provides an auxiliary USB power source for charging devices and optional LED lights for improved visibility
- **RAM® Mounts Self-Leveling Cup Holder** uses gimbal system to keep drinks level during travel or when stopped on inclines
- **Active Controls ReCharge™ with Socket** eliminates the challenges of connecting and disconnecting the battery charger by providing an easy-to-reach adapter

Quickie Pulse Specifications

| SPECIFICATIONS | PULSE 5 BC | PULSE 6 BC | PULSE 5 CC | PULSE 6 CC | PULSE 6 SC | PULSE 6 MPC |
|-----------------------|--------------------|--------------------|-----------------------|-----------------------|-----------------------------------|------------------|
| Drive Wheel Position | Mid | Mid | Mid | Mid | Mid | Mid |
| HCPCS Code | K0848 | K0848 | K0849 | K0849 | K0856 | K0861 |
| Speed | 5.0 MPH | 6.0 MPH | 5.0 MPH | 6.0 MPH | 6.0 MPH | 6.0 MPH |
| Overall Weight | 235 lbs. | 235 lbs. | 250 lbs. | 250 lbs. | 260 lbs. | 345 lbs. |
| Overall Width | 24" | 24" | 24" | 24" | 24" | 24" |
| Overall Length | 35" | 35" | 35" | 35" | 35" | 35" |
| User Weight Capacity | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. |
| Seat Widths | 12" to 20" | 12" to 22" | 16", 18", 20" | 16", 18", 20" | 12" to 22" | 14" to 22" |
| Seat Depths | 12" to 20" | 12" to 22" | 15" to 22" | 15" to 22" | 12" to 22" | 16" to 20" |
| Seat To Floor Height | 16" to 20" | 16" to 20" | 16" to 18" (from pan) | 16" to 18" (from pan) | 16.3" to 21.6" | 16.3" to 18.3" |
| Seating Frame Options | ASAP II Rehab Seat | ASAP II Rehab Seat | Captain's Seat | Captain's Seat | ASAP II Rehab Seat or 3.7 Recline | 3.7 Recline |
| Transit Approved | ISO | ISO | WC-19* | WC-19* | Tilt: ISO/Recline: WC-19* | WC-19* |
| Warranty | Limited/Lifetime | Limited/Lifetime | Limited/Lifetime | Limited/Lifetime | Limited/Lifetime | Limited/Lifetime |

* Conforms with ANSI/RESNA WC-19 standards. Please refer to product Owner's Manual for additional information. Owner's Manuals are available for download at www.sunrisemedical.com

Live Without Limits

The Quickie S-636 features innovative, Group 3 performance and versatile seating and positioning options to accommodate both general and specialized rehabilitation needs.

QUICKIE® S-6 Series



Available in Six Group 3 and 4 Models:


- Quickie S-636 BC (K0848): No Power Option S-636 base with Rehab seat
- Quickie S-636 SC (K0856): Single Power Option (tilt or recline) S-636 base with Rehab or Recline seat
- Quickie S-636 MPC (K0861): Multi Power Option (tilt and recline) S-636 base with Rehab Recline seat
- Quickie S-646 SE BD (K0868): No Power Option S-646 SE base with Rehab seat
- Quickie S-646 SE SD (K0877): Single Power Option (tilt or recline) S-646 SE base with Rehab or Recline seat
- Quickie S-646 SE MPD (K0884): Multi Power Option (tilt and recline) S-646 SE base with Rehab Recline seat



Hybrid Base with Dynamic Suspension

- Linix 6.5 mph or 8.5 mph, 4-pole standard motor package for all models
- Short wheel base with hybrid wheel placement for mid-wheel maneuverability and rear-wheel all-terrain performance
- Adjustable dampening, independent wheel suspension that provides a smooth ride, even over obstacles
- PosiDrive™ technology increases the overall strength of the drive train
- Optional LED headlight, blinker and tail light package for safety
- Painted frame available in 23 different color options

Seating Options



Rehab Seat

- TRAX System with easy dovetail track adjustments of the seat and back
- Extensive seating and positioning options
- Articulating vent tray option

Power Recline Seat

- Wide range of seating and positioning options
- Low seat to floor options with and without Power CG Tilt
- Headrest design is easily accessed and does not require tools for removal



Power Options

- 53° of Power CG Tilt for optimized positioning
- Power Recline with Power Shear Reduction
- Power Elevating Legrests with Articulation

Quickie S-6 Series Specifications

| SPECIFICATIONS | S-636 BC | S-636 SC | S-636 MPC | S-646 SE BD | S-646 SE SD | S-646 SE MPD |
|----------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|
| HCPCS Code | K0848 | K0856 | K0861 | K0868 | K0877 | K0884 |
| Drive Wheel Position | RWD | RWD | RWD | RWD | RWD | RWD |
| Speed | 6.5 MPH* | 6.5 MPH* | 6.5 MPH** | 8.5 MPH | 8.5 MPH | 8.5 MPH |
| Overall Weight | 128 lbs. w/o batteries |
| Overall Width | 27.75" | 27.75" | 27.75" | 25.75" | 25.75" | 25.75" |
| Overall Length | 33.5" | 33.5" | 33.5" | 33.5" | 33.5" | 33.5" |
| User Weight Capacity | 300 lbs. |
| Seat Widths | 14" to 24" |
| Seat Depths | 14" to 22" |
| Front Seat To Floor Height | 18.5" to 20.5" | 18.5" to 20.5" | 18.5" to 20.5" | 18" to 21" | 18.5" to 20.5" | 18.5" to 20.5" |
| Seating Frame Options | Rehab | Rehab or Rehab Recline | Rehab Recline | Rehab | Rehab or Rehab Recline | Rehab Recline |
| Power Seating Options | N/A | Til or Recline | Tilt and Recline | N/A | Tilt or Recline | Tilt and Recline |

* Optional 5 mph and 8.5 mph motor packages also available

** Optional 8.5 mph motor package also available

Live Without Limits

The powerful and compact Quickie P-222 SE is fast, simple, and extremely maneuverable indoors and out and features a quick-release power pack and optional manual conversion kit for convenience and transportability.

QUICKIE® P-222 SE



Available in Two Group 4 Models:

- Quickie P-222 SE BD (K0868): No power option P-222 SE base with Rehab seat
- Quickie P-222 SE SD (K0877): Single power option (tilt) P-222 SE base with Rehab seat



Fast and Maneuverable Base

- Linix 8.5 mph, 4-pole standard motor package
- PosiDrive™ technology increases the overall strength of the drive train and features easy-to-service hubs
- 26" turning radius for excellent maneuverability
- 4" curb climbing capability



Seating and Power Options

Rehab Seat

- Simple seating and positioning options mix
- Swing-away and fixed hanger options
- Optional 53° of CG tilt to optimize positioning



Style Options

- Optional manual wheelchair conversion kit for user self propulsion
- Optional LED light package includes headlights, blinkers and tail light for safety
- Painted frame available in 23 different color options

Quickie P-222 SE Specifications

| SPECIFICATIONS | P-222 SE BD | P-222 SE SD |
|----------------------------|------------------------|------------------------|
| HCPCS Code | K0868 | K0877 |
| Drive Wheel Position | RWD | RWD |
| Speed | 8.5 MPH | 8.5 MPH |
| Overall Weight | 115 lbs. w/o batteries | 115 lbs. w/o batteries |
| Overall Width | 24.5" | 24.5" |
| Overall Length | 31.5" w/o footplates | 31.5" w/o footplates |
| User Weight Capacity | 350 lbs. | 350 lbs. |
| Seat Widths | 14" to 24" | 14" to 24" |
| Seat Depths | 14" to 20" | 14" to 20" |
| Front Seat To Floor Height | 18" to 20" | 17.5" to 20.5" |
| Seating Frame Options | Rehab | Rehab |
| Power Seating Options | N/A | Tilt |

PG Drives Technology

Option availability varies by chair and model. Please see order form for detailed configuration information.



Switch-It Head Array

VR2 Non-Expandable Controls

- Industry standard non-expandable electronics
- Dependable, highly flexible, and simple to use
- 4-way toggle switch for operation of tilt and power legrests without the need for expandable controls
- Optional Attendant Control for operation by caregiver

R-net Expandable and Non-Expandable Controls

- 8 uniquely programmable profiles for easy transitions between environments
- Optional assignable buttons for convenient control of common functions
- Versatile joystick options
- Optional LED joystick
- Optional joystick w/color display
 - Color Screen for clear operation of functions
 - On-board programming capabilities
 - Available with mono and stereo jacks for adding external actuator and ON/OFF switch control
 - Available with Toggle, Speed Potentiometer, and stereo jacks for quad-friendly controls
- Expandable OMNI Input Control Module
 - Infra Red (IR) functions for operation of IR controlled devices such as TVs and DVD players
 - Bluetooth Mouse-Mover for PC mouse control
 - On-board programming capabilities
- Full specialty control and switch input capabilities
- Optional Attendant Control for operation by caregiver

Back Upholstery Options



QUICKIE® Tension Adjustable Upholstery

- Available with Ballistic Nylon or 3DX Microclimatic fabric
- 1" width adjustment




SPOT Curved Back

- Solid Support
- Optional preassembled Whitmyer headrests
- Optional AES positioning components



JAY® J3™ Backrest

- Two-point or four-point quick-release hardware
- Range of back heights and contours to meet user needs
- Can mount Whitmyer headrests or lateral positioning supports



SPO Backrest

- JAY or ART SPO Backrests are designed for the unique needs of a client in power tilt.
- Provides postural management in a package that can be customized in the field.
- Shape, height, width and thickness of the foam customizable by user.
- ART SPO Back offers custom shapes: i-back, scapular cut-out and custom rectangular.

Electronics

Option availability varies by chair and model. Please see order form for detailed configuration information.

Drive Controls



VR2 Non-Expandable Controller



R-net Expandable Controller



R-net 120amp Non-Expandable/Expandable Controller



R-net Environmental Control Output Module (ECM)



R-net Bluetooth® Mouseover

Joysticks



Non-expandable VR-2 Joystick



R-net Expandable LED Joystick



Joystick w/Stereo Jacks & Color Display



Joystick w/Toggle, Speed Potentiometer, Stereo Jacks & Color Display



R-net OMNI Input Control Module

Joystick Mounts



Fixed



Swing-Away



VR-2



R-net

Attendant Controls

Joystick Handles



Ball



Foam Ball



T-Handle



Goalpost

Option availability varies by chair and model. Please see order form for detailed configuration information.

Switch Inputs for Joysticks or Omni



Ribbon Switch



Buddy Button Switch



Disc Switch



Wobble Switch



Micro Light Switch



Egg Switch - Black



Stereo to 2 Mono Splitter




Switch-It® Mini 1 Button on Mono Plug



Switch-It® Mini 2 Button on 2 Mono Plugs



Switch-It® Mini 2 Button on 1 Stereo Plug



Switch-It® Mini 4 Button on DB9 Plug

Proportional Controls - Chin Control



Switch-It® Micro-Pilot Bib Chin Control Pkg, Standard 6" or Large 8"



Switch-It® Micro-Guide Bib Chin Control Pkg, Standard 6" or Large 8"



Quickie® Proportional Chin Control Package



ASL® Proportional Micro Extremity Joystick Chin/Lip Control Package

Proportional Controls - Hand Control



ASL® Proportional Micro Extremity Joystick Midline with Gatlin Package



ASL® Mushroom Joystick Midline with Gatlin Package



ASL® Compact Joystick with Integrated Buttons Midline with Gatlin Package



ASL® Micro Mini Joystick with Hand Pad Midline with Gatlin Package



Switch-It® Micro-Pilot Joystick with Bullet Tray Package



Switch-It® Micro-Guide Joystick with Bullet Tray Package

Switch Control



ASL® Proximity Head Array Package 3-Switch (Adult or Pediatric)



Switch-It® Head Array Package 3-Switch Curved (Long or Stubby)



Switch-It® Head Array 5-Switch Flat



Quickie® Sip 'N' Puff System


Multiple Switches



Switch-It® 3 Naked Proximity Switches on DB9 Connector



Switch-It® 4 Naked Proximity Switches on DB9 Connector



Switch-It® 5 Naked Proximity Switches on DB9 Connector



ASL® 5-Switch Adaptor Box on DB9 Connector

Programmers



VR-2 Hand Held Programmer



R-net PC Programmer



DTT VR2 and R-net

Base Options


Option availability varies by chair and model. Please see order form for detailed configuration information.



Knobby Drive Wheels



4-Point Transit Tie-Downs



Fixed Vent Tray



Gp 24 Batteries with Std Posts



Gp 24 Batteries with Flat Tops



22NF Batteries with Std Posts



22NF Batteries with Flat Tops



Gp 34 Batteries with Flat Tops



8A Dual Mode (Off-Board) CTE

Accessories

Option availability varies by chair and model. Please see order form for detailed configuration information.



RAM® Cup Holder



RAM® Phone Mount



RAM® Tablet Mount



Active Controls Power Buddy



Active Controls ReCharge



Switch-It® Remote Stop Switch



Auto Style Buckle



Belt 2" Aircraft Buckle



Belt 2" Aircraft Buckle Padded



O2 Holder



Crutch Holder



Walker Holder



Tuffak Tray Table



Black Utility Backpack



Black/Charcoal Kid's Backpack

Armrests

Option availability varies by chair and model. Please see order form for detailed configuration information.



**Dual Post Flip Back
Height Adjustable**



Flip Back



**Standard Cantilever with 2 Point
Locking Post**



Standard Cantilever



Low Cantilever



Standard Single Post



Reclining Standard Height



Cantilever Reclining

Armpads



Waterfall Armpad



Standard Armpad



Reclining Armpad



Ottobock® Reclining Armpad



Ottobock® Flat Hand Pad

Positioning Options



***Whitmyer® Plush
with Lynx hardware**



***Whitmyer® Plush
with Cobra hardware**



***Whitmyer® Plush
with Axis hardware**



Chest Strap



Lateral Thoracic Supports



JAY® J3™ Adjustable Lateral



Swing-Away Knee Adductor



**Long Thigh Support
with mesh pocket**

**Various Whitmyer head supports are available with the hardware above.*

Legrests

Option availability varies by chair and model. Please see order form for detailed configuration information.



65° & 70° Swing-Away



75° Tapered Swing-Away



90° Swing-Away



Manual ELR Swing-Away



Fixed Center Mount



Centermount Power ELR with Articulation



Swing Away Power ELR with Articulation



Kids Angle Adjust with Leg Strap

Footplates



Composite with Heel Loop



Adult Angle Adjust with Heel Loop



Kids Angle Adjust with Leg Strap



Aluminum



One Piece Angle Adjust Footplate



Pediatric Footplate



Centermount Single/Dual Footplate



www.built4me.us

Product Customizations

Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint!

Featured Power Chair Customizations

- 24" Seat width on QM-715 HD
- Captain's Seat on QM-710

JAY[®] Seating Selection

JAY CUSHIONS



Basic
E2601, E2602



BasicPRO™
E2601, E2602



GO
E2603, E2604



Soft Combi™ P
E2605, E2606



Care
E2607, E2608



Zip™
E2607



Ion™
E2607, E2608



Easy
E2607, E2608



Union™
E2307, E2608



X2™
E2622



J2™
E2622, E2623,
E2624, E2625



Fusion™
E2622, E2623



J3™
E2622, E2623,
E2624, E2625

JAY BACKS



Basic
E2611



GO
E2611, E2612



Care
E2615



Focus Point
E2620



Zip™
E2613



J3™ Posterior
E2613, E2614



**J3™ Posterior
& Lateral**
E2615



**J3™ Posterior
& Deep Lat/Cont**
E2620